



Session Supplies for Open Colonic Devices

Email for New Account Application



Cleansing Nozzle: Curved Bulb Tip - RX


Flexible Tubing FDA Food Grade

No Leaching of dangerous polymers, chemicals or resins into human body.

50 + Ft. Roll
(You cut to size)



Plastic Handle Mop
Mop-on-rope (snake)



Petrolatum Water Base Jelly
Package Colors may Change

Degreaser Concentrated

Cleans View Tube

- Gallon • LS-36

“CDC”. Disinfectant
FDA Approved
(Neutralizes HIV!)
Concentrated
LS-34 • Gallon



Maintenance Parts



Carbon Filter

Hot Water Use.
MUST BE REPLACED MONTHLY TO PREVENT BACTERIA

Rated: 5 Micron
Reduces: Chlorine, Odor, Sediment



UV Bulb - Replaced Annually

Available for LIBBE's Built After 2002

CONNECTOR NOZZLE

Flex Tubing connection can wear out - Barbed End
Replace when needed!




Tummy Warmer

Reusable over & over
Simply boil after using.
Stays warm 60 minutes



Gloves Nitrile

200 PER BOX
Med, Large, XLarge





EZWAX

For waxing the LIBBE
Fiberglass Base / will look new!



Wedge



Cushion Set (2 Pc.)

Tiller MIND BODY, Inc.
10911 West Avenue
San Antonio, Texas USA
78213
Office: 210 308.8888
Email: info@colonic.net



Made in USA